

TT2 ANESTHESIA SYSTEM

**Table Top 2* Anesthesia System
with built in fresh gas dual delivery valve.**

FEATURES

- ✓ Compact 12 inch footprint.
- ✓ Compatible with all vaporizer mounts.
- ✓ Designed for ease of use, setup and portability.
- ✓ Flowmeter with mechanical stop flow control knob, 0-4 LPM.
- ✓ Customizable with the ability to add additional accessories.


QUALITY

- ✓ Machined from anodized aluminum for durability and chemical resistance.
- ✓ Ability to utilize one or both fresh gas delivery outlets without disconnecting supply lines.
- ✓ 5 year warranty.

INCLUDES

- ✓ Vaporizer inlet and outlet endcaps.
- ✓ Color coded fresh gas delivery lines with common outlet and 15 mm adaptor ends (2).

**Vaporizer and other accessories not included.*



SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625 ✉ info@somniscientific.com

🌐 www.somniscientific.com



*Proudly machined and
assembled in the USA*