



INDUCTION CHAMBERS

ACTIVE VACUUM AND PASSIVE STYLES

The SOMNI line of induction chambers are designed with slide top lids that allow for one handed operation. Unlike a “hinge lid” the slide top keeps gases from being “drawn” into the workplace when the chamber is opened.



QUALITY

- ✓ Constructed of durable, shatter resistant polycarbonate.
- ✓ Can be safely cleaned with common clinical disinfectants.

SAFETY

- ✓ Prevents unnecessary exposure to personnel and environment (NIOSH 2ppm) when used with active waste gas removal system EPS-3 or W.A.G. filter (active chamber only).

CONVENIENT

- ✓ One handed operation.
- ✓ Color coded inlet and outlet ports.
 - Compatible with 19mm tubing & 15mm adapters.

SPECIFICATIONS

ACTIVE

Large	11.5 in L	6 in W	6 in H
-------	-----------	--------	--------

PASSIVE

Small	4.25 in L	4.25 in W	4.25 in H
Medium	9.25 in L	5.25 in W	5.25 in H
Large	11.5 in L	6 in W	6 in H

SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625 ✉ info@somniscientific.com

🌐 www.somniscientific.com



Proudly machined and assembled in the USA