

# SOMNI W.A.G. WASTE ANESTHETIC GAS COLLECTION FILTER

The SOMNI W.A.G. Waste Anesthetic Gas Collection Filter is an activated charcoal filter developed to absorb 200 grams of halogenated agent. The W.A.G. filter can be incorporated into an inhalant anesthetic system to passively and safely remove waste anesthetic gas at higher flowrates and from multiple sources.



## QUALITY

- ✓ Absorbs up to 200 grams of halogenated agent.
- ✓ Inlet dimensions (15mm ID, 19/22mm OD) accommodate multiple tubing options.

## SAFETY


- ✓ Prevents unnecessary exposure to personnel and environment (NIOSH 2ppm) when replaced regularly.
- ✓ Use where evacuation to the outside is not an option.
- ✓ Filter can be placed on a flat, stable surface during use.

## CONVENIENT

- ✓ Easily adapted for use with multiple sources of WAG .
- ✓ Disposable.

## SPECIFICATIONS

Weight	Approximate 2 lbs. / 1 kg
Height	9 in / 18 cm
Diameter	5 in / 23 cm
Count	4 Pack



**SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.**

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625    ✉ info@somniscientific.com

🌐 www.somniscientific.com



*Proudly machined, assembled  
and calibrated in the USA*