



AIRSEP® VISIONAIRE O2 CONCENTRATOR

QUALITY

- ✓ Cabinet is made of Super High Impact Polystyrene.
- ✓ Brass gas outlet and anodized black aluminum gas nozzle.
- ✓ Durable, maintenance free design.
- ✓ Consistent supply of high concentration supplemental oxygen.

ECONOMICAL

- ✓ Low total cost of ownership.
- ✓ Low consumption power of 290 W at all flow rates.
- ✓ Eliminates the requirement of delivered oxygen.

SAFE

- ✓ Lockable flowmeter.
- ✓ Safety pressure relief valve calibrated to 280 kPa (40psig).
- ✓ Power failure, no flow, high and low pressure alarms.
- ✓ LED Oxygen monitor visual alarm when concentration drops below 82%.

CONVENIENT

- ✓ Top and Side handles.
- ✓ Casters for portability.
- ✓ Lightweight.
- ✓ Smooth quiet operation.
- ✓ Optional humidifier bottle.

SPECIFICATIONS

Weight	30 lbs. / 13.6 kg
Height	20.8 in / 52.8 cm
Depth	11.5 in / 29.2 cm
Width	14.1 in / 35.8 cm
Flow Meter Capacity	1-5 LPM (+/- 10% of indicated setting or 200mls, whichever is greater)
Oxygen Concentration	90% (+ 5.5% / -3%)
Sound	40 - 45 dBa

SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625 ✉ info@somniscientific.com

🌐 www.somniscientific.com